

# SkinRevive(LITE)

Boosts Collagen Naturally

the beauty specialists



# SkinRevive(LITE)

SKINREVIVE™ uses photomodulation with a 633nm red light source to naturally rejuvenate the skin and counter the signs of photoageing.

The targeted red light source boosts natural collagen production in the skin to help improve skin texture and elasticity.

This results in a firmer, smoother and more youthful appearance.





# SkinRevive(LITE)

## Features

- Uses photomodulation with a 633nm red light

## Benefits

- Reduces Fine Lines and Wrinkles
- Firms and Tones
- Promotes Collagen Production
- Rejuvenates Skin Tone





# SkinRevive(LITE)

Thoroughly cleanse your face before treatment.

## **STEP 1**

Briefly press and hold the ON/OFF button to switch on.

## **STEP 2**

Gently apply the treatment head to the skin.

## **STEP 3**

Gently move the treatment head over the surface of the skin focusing on areas showing signs of ageing. As a guide, each half of the face should be treated for approximately 10 minutes.

## **STEP 4**

To switch off, press and hold the ON/OFF button. Wipe the treatment head clean with an antiseptic wipe ready for your next treatment.

You can follow with your normal beauty routine immediately after treatment.

# SkinRevive(LITE)

## Contents:

- Skin Revive unit
- Charging Dock
- Power Adapter

- **Program Time:** 20 mins
- **Power Levels:** 1
- **Output:** 633nm
- **Packed Product Dimensions:** (height x width x depth) 225 x 130 x 85 mm
- **Un-packed Product Dimensions:** (height x width x depth) 200 x 58 x 58 mm
- **AC/DC Adaptor:** No
- **Battery Operated/Type:** Rechargeable
- **Manufacturer's Age Recommendation:** 14 and over
- **Packed Weight (Kg):** 0.436
- **Product Net Weight (Kg):** 0.264
- **Manufacturer's Guarantee:** 2 years
- **Barcode:** 5019487083531

# Frequently Asked Questions About **SkinRevive** (LITE)

## **How quickly will I get results?**

Clinical trials of photomodulation have shown improvements after just 8 treatments over 4 weeks for people with mild to moderate photoageing.

## **Can I use whilst pregnant?**

Yes.

## **Is it safe for use around the eyes?**

Yes, but be careful to avoid looking directly at the LEDs. If you glance accidentally at the LEDs, no harm will occur. We recommend keeping your eyes closed during treatment to reduce the effects of glare. If you suffer with any eye disorders then please consult your specialist before use.

## **Does it work for all skin types?**

Yes.

## **Can I use it on other areas of the body?**

Because this is a natural rejuvenation process, you can also use it on your neck, hands or areas of fatigued skin.

## **Can I use it with creams?**

The SKINREVIVE is intended to be a natural treatment - no creams are necessary.